

1.0 L and 1.5 L Prefilled Closed Enteral Feeding Containers May Be Hung for Up to 48 Hours

1.0 or 1.5 Liter Containers



+

Closed Feeding Spike System



Procedure:

- **Continuous/ Uninterrupted Infusions – Max hang time – 48 hours**
 - Spike bag once and use clean technique.
- **Interrupted infusions – Max hang time – 24 hours**
 - Interruption occurs if connector feeding set is disconnected from the patient or the container is spiked more than once.

Benefits:

- Saves nursing time.
- Minimizes risk of touch contamination.
- Conserves connector feeding sets which are in short supply nationally due to COVID-19.