

# Visitor Mask Assembly Instructions

1. For providing to hospital visitors
2. Typical assembly ~60 seconds or less
3. Patience during assembly minimizes damages
4. Choose a clean area for assembly
5. Wear mask and gloves during assembly

# Mask Unassembled



MedStar Health

Knowledge and Compassion **Focused on You**

# Step One:

- **Fold mask in half**
- **Place bottom center tab into top center slit**



# Step Two:

---

- Thread tie from one side through three slits on same side
- Pull together through slits





# Step Three:

- Repeat Step Two for the other side
- Thread tie from other side through three slits on same side
- Pull together through slits



# Step Four:

---

- **Mask is ready for use**
- **Bottom ties are tied around the neck**
- **Top ties are tied around the head, ensuring that nose and mouth are fully covered**

